

# PolyFlush

## Pre-filled NaCl syringe for flushing

PolyFlush pre-filled syringes are administered before and after starting IV medication as per standard protocols to ensure that IV catheters are clear and free from obstruction. Maintain patency within the catheter so that blood is prevented from backfilling into the catheter and clotting.

**Avoids catheter damage:** Polyflush prefilled syringe offers all the syringes in consistent 10 ml syringe diameters but adapted length according to the content in ml. Smaller diameter syringes generate greater amounts of pressure than larger diameter syringes and can risk catheter damage.

### Features:

- Contains 0.9% (9 mg/ml) of isotonic sodium chloride (NaCl).
- Ready to use, limits number of steps, fewer connections, and less packaging waste.
- Polyflush is a pre-labelled prefilled saline syringe which may reduce the risk of medication error.
- PolyFlush prefilled syringes minimizes cross contamination.
- The Polyflush pre-filled syringes are available in two different versions, Single Sterile (SS) or Double Sterile (DS)\*.
- Available in sizes of 3ml, 5ml and 10ml.
- Gamma sterilized: reduces the risk of smell and taste sensation.
- PVC free, DEHP free & Latex free

\*When inserting intravascular catheters/sterile dressing, a Double Sterile prefilled syringe should be used.

**Single Sterile, SS**  
3ml 30 pcs/box  
REF: 90311



**Double Sterile, DS**  
3ml 30 pcs/box  
REF: 90317



**Single Sterile, SS**  
5ml 30 pcs/box  
REF: 90319




**Double Sterile, DS**  
5ml 30 pcs/box  
REF: 90316



**Single Sterile, SS**  
10ml 30 pcs/box  
REF: 90318



**Double Sterile, DS**  
10ml 30 pcs/box  
REF: 90315



PolyFlush is a medical device only intended for flushing intravenous catheters. Therefore, not to be used for mixing substances, diluting medicine, or used as intravenous treatment.